

# TAVA SURGICAL<sup>®</sup>

## HIGH SPEED BURS



PRECISION. PERFORMANCE.



Your Value Provider



**TAVA SURGICAL<sup>®</sup>**  
SURGICAL POWER & ACCESSORIES




TAVA Surgical's high Speed burs are offered at significant cost savings and are fully compatible with the original manufacturer's equipment. Our high speed burs are designed and manufactured to very rigid specifications and are subjected to extensive engineering analysis and strenuous testing programs to ensure exceptional product performance.

## TAVA Surgical® High Speed Burs Compatible with Medtronic® Midas Rex® Legend® Systems

TAVA Surgical® Order Number	Replaces	Description	Head Diameter (mm)	Head Length (mm)	TAVA Surgical® High Speed Burs Compatible with the following Attachments
					
<b>BR9AC60</b>	<i>9AC60</i>	Acorn - Fluted	6.0	7.7	AA09S, AS09, AS09S, AVA09, AVS09
<b>BR9AC75</b>	<i>9AC75</i>	Acorn - Fluted	7.5	9.4	AA09S, AS09, AS09S, AVA09, AVS09
<b>BR9AC90</b>	<i>9AC90</i>	Acorn - Fluted	9.0	11.4	AA09S, AS09, AS09S, AVA09, AVS09
<b>BR14AC60</b>	<i>14AC60</i>	Acorn - Fluted	6.0	7.7	AA14, AA14S, AS14, AS14S, AVA14, AVS14
<b>BR14AC75</b>	<i>14AC75</i>	Acorn - Fluted	7.5	9.4	AA14, AA14S, AS14, AS14S, AVA14, AVS14
					
<b>BR7BA20</b>	<i>7BA20</i>	Ball - Fluted	2.0	N/A	AA07, AS07 AVA07, AVS07
<b>BR7BA50</b>	<i>7BA50</i>	Ball - Fluted	5.0	N/A	AA07, AS07 AVA07, AVS07
<b>BR9BA30</b>	<i>9BA30</i>	Ball - Fluted	3.0	N/A	AA09S, AS09, AS09S, AVA09, AVS09
<b>BR9BA40</b>	<i>9BA40</i>	Ball - Fluted	4.0	N/A	AA09S, AS09, AS09S, AVA09, AVS09
<b>BR9BA50</b>	<i>9BA50</i>	Ball - Fluted	5.0	N/A	AA09S, AS09, AS09S, AVA09, AVS09
<b>BR9BA60</b>	<i>9BA60</i>	Ball - Fluted	6.0	N/A	AA09S, AS09, AS09S, AVA09, AVS09
<b>BR9BA75</b>	<i>9BA75</i>	Ball - Fluted	7.5	N/A	AA09S, AS09, AS09S, AVA09, AVS09
<b>BR10BA30</b>	<i>10BA30</i>	Ball - Fluted	3.0	N/A	AA10, AA10S, AS10, AS10S, AVA10, AVS10
<b>BR10BA40</b>	<i>10BA40</i>	Ball - Fluted	4.0	N/A	AA10, AA10S, AS10, AS10S, AVA10, AVS10
<b>BR10BA50</b>	<i>10BA50</i>	Ball - Fluted	5.0	N/A	AA10, AA10S, AS10, AS10S, AVA10, AVS10
<b>BR14BA30</b>	<i>14BA30</i>	Ball - Fluted	3.0	N/A	AA14, AA14S, AS14, AS14S, AVA14, AVS14
<b>BR14BA40</b>	<i>14BA40</i>	Ball - Fluted	4.0	N/A	AA14, AA14S, AS14, AS14S, AVA14, AVS14
<b>BR14BA50</b>	<i>14BA50</i>	Ball - Fluted	5.0	N/A	AA14, AA14S, AS14, AS14S, AVA14, AVS14






## TAVA Surgical® High Speed Burs Compatible with Medtronic® Midas Rex® Legend® Systems

TAVA Surgical® Order Number	Replaces	Description	Head Diameter (mm)	Head Length (mm)	TAVA Surgical® High Speed Burs Compatible with the following Attachments
					
<b>BR9BA50D</b>	<i>9BA50D</i>	Ball - Diamond	5.0	N/A	AA09S, AS09, AS09S, AVA09, AVS09
<b>BR10BA50D</b>	<i>10BA50D</i>	Ball - Diamond	5.0	N/A	AA10, AA10S, AS10, AS10S, AVA10, AVS10
<b>BR14BA50D</b>	<i>14BA50D</i>	Ball - Diamond	5.0	N/A	AA14, AA14S, AS14, AS14S, AVA14, AVS14
<b>BR9BA30DC</b>	<i>9BA30DC</i>	Ball - Diamond Coarse	3.0	N/A	AA09S, AS09, AS09S, AVA09, AVS09
<b>BR9BA60DC</b>	<i>9BA60DC</i>	Ball - Diamond Coarse	6.0	N/A	AA09S, AS09, AS09S, AVA09, AVS09
<b>BR14BA60DC</b>	<i>14BA60DC</i>	Ball - Diamond Coarse	6.0	N/A	AA14, AA14S, AS14, AS14S, AVA14, AVS14
<b>BR9BA60DX</b>	<i>9BA60DX</i>	Ball - Diamond Extra Coarse	6.0	N/A	AA09S, AS09, AS09S, AVA09, AVS09
<b>BR10BA40DX</b>	<i>10BA40DX</i>	Ball - Diamond Extra Coarse	4.0	N/A	AA10, AA10S, AS10, AS10S, AVA10, AVS10
					
<b>BR9CY50</b>	<i>9CY50</i>	Cylinder - Fluted	5.0	7.9	AA09S, AS09, AS09S, AVA09, AVS09
<b>BR10CY40</b>	<i>10CY40</i>	Cylinder - Fluted	4.0	6.2	AA10, AA10S, AS10, AS10S, AVA10, AVS10
<b>BR14CY50</b>	<i>14CY50</i>	Cylinder - Fluted	5.0	7.9	AA14, AA14S, AS14, AS14S, AVA14, AVS14
<b>BR14CY60</b>	<i>14CY60</i>	Cylinder - Fluted	6.0	7.9	AA14, AA14S, AS14, AS14S, AVA14, AVS14
					
<b>BR14CY65DX</b>	<i>14CY65DX</i>	Cylinder - Diamond Extra Coarse	6.5	9.8	AA14, AA14S, AS14, AS14S, AVA14, AVS14



## TAVA Surgical High Speed Burs Compatible with Medtronic® Midas Rex® Legend® Systems

TAVA Surgical® Order Number	Replaces	Description	Head Diameter (mm)	Head Length (mm)	TAVA Surgical® High Speed Burs Compatible with the following Attachments
					
<b>BR8MH17</b>	<i>8MH17</i>	Match Head - Fluted	1.7	2.4	AS08, AVS08
<b>BR9MH30</b>	<i>9MH30</i>	Match Head - Fluted	3.0	3.8	AA09S, AS09, AS09S, AVA09, AVS09
<b>BR10MH22</b>	<i>10MH22</i>	Match Head - Fluted	2.2	3.8	AA10, AA10S, AS10, AS10S, AVA10, AVS10
<b>BR10MH30</b>	<i>10MH30</i>	Match Head - Fluted	3.0	3.8	AA10, AA10S, AS10, AS10S, AVA10, AVS10
<b>BR10MH30-MN</b>	<i>10MH30-MN</i>	Match Head - Fluted (MedNEXT® Type)	3.0	3.8	AA10, AA10S, AS10, AS10S, AVA10, AVS10
<b>BR14MH30</b>	<i>14MH30</i>	Match Head - Fluted	3.0	3.8	AA14, AA14S, AS14, AS14S, AVA14, AVS14
<b>BR15MH22</b>	<i>15MH22</i>	Match Head - Fluted	2.2	3.8	AA15, AS15 AVA15, AVS15
					
<b>BR9MH30D</b>	<i>9MH30D</i>	Match Head - Diamond	3.0	3.8	AA09S, AS09, AS09S, AVA09, AVS09
<b>BR10MH30D-MN</b>	<i>10MH30D-MN</i>	Match Head - Diamond (MedNEXT® Type)	3.0	4.1	AA10, AA10S, AS10, AS10S, AVA10, AVS10
<b>BR14MH30D</b>	<i>14MH30D</i>	Match Head - Diamond	3.0	3.8	AA14, AA14S, AS14, AS14S, AVA14, AVS14
					
<b>BR10OV40</b>	<i>10OV40</i>	Oval	4.0	6.4	AA10, AA10S, AS10, AS10S, AVA10, AVS10
<b>BR14OV40</b>	<i>14OV40</i>	Oval	4.0	6.4	AA14, AA14S, AS14, AS14S, AVA14, AVS14



## TAVA Surgical® High Speed Burs Compatible with Medtronic® Midas Rex® Legend® Systems

TAVA Surgical® Order Number	Replaces	Description	Head Diameter (mm)	Head Length (mm)	TAVA Surgical® High Speed Burs Compatible with the following Attachments
<b>BR8TA11</b>	<i>8TA11</i>	Tapered	1.1	6.4	AS08, AVS08
<b>BR14TA30</b>	<i>14TA30</i>	Tapered	3.0	35.1	AA14, AA14S, AS14, AS14S, AVA14, AVS14
<b>BRF1/8TA15</b>	<i>F1/8TA15</i>	Tapered	1.5	11.1	AF01, AF01R, AS08, AVS08
<b>BRF2/8TA23</b>	<i>F2/8TA23</i>	Tapered	2.3	15.9	AF02, AF02R, AS08, AVS08
<b>BR8TD114</b>	<i>8TD114</i>	Twist Drill	1.1	4.0 Stop Depth	AS08, AVS08
<b>BR8TD116</b>	<i>8TD116</i>	Twist Drill	1.1	6.0 Stop Depth	AS08, AVS08
<b>BR8TD158</b>	<i>8TD158</i>	Twist Drill	1.5	8.0 Stop Depth	AS08, AVS08



## TAVA Surgical® High Speed Burs Compatible with Medtronic® Midas Rex® Legend® Systems


TAVA Surgical® Order Number	Replaces	Description	Head Diameter (mm)	Head Length (mm)	TAVA Surgical® High Speed Burs Compatible with the following Attachments
<b>BRMC254</b>	<i>MC254</i>	Metal Cutter	25.4	1.0	ASMC
<b>BRMC30</b>	<i>MC30</i>	Metal Cutter	3.0	18.3	ASMC

TAVA Surgical's high speed burs are offered at significant cost savings and are fully compatible with the original manufacturer's equipment. Our high speed burs are designed and manufactured to very rigid specifications and are subjected to extensive engineering analysis and strenuous testing programs to ensure exceptional product performance.

## TAVA Surgical® High Speed Replacement Burs





### Compatible with Stryker® TPS™, Medtronic® Midas Rex® Classic I & II and Linvatec® UltraPower™ Systems



Order Number	Description	Unit	Overall Length	Head Diameter	Head Length
					
<b>KM-M330U</b>	Round Fluted	mm	83.0	3.0	3.0
		in	3.27	0.12	0.12
<b>KM-M340U</b>	Round Fluted	mm	84.0	4.0	4.0
		in	3.31	0.16	0.16
<b>KM-M35U</b>	Round Fluted	mm	86.0	5.1	5.1
		in	3.39	0.20	0.20
<b>KM-M32U</b>	Round Fluted	mm	89.0	7.6	7.6
		in	3.50	0.30	0.30
<b>KM-M31U</b>	Round Fluted	mm	91.0	9.1	9.1
		in	3.58	0.36	0.36
<b>KM-AM330U</b>	Round Fluted	mm	137.0	3.0	3.0
		in	5.39	0.12	0.12
<b>KM-AM340U</b>	Round Fluted	mm	137.0	4.0	4.0
		in	5.39	0.16	0.16
<b>KM-AM35U</b>	Round Fluted	mm	140.0	5.1	5.1
		in	5.51	0.20	0.20
<b>KM-AM33U</b>	Round Fluted	mm	140.0	6.0	6.0
		in	5.51	0.24	0.24
<b>KM-AM32U</b>	Round Fluted	mm	143.0	7.6	7.6
		in	5.63	0.30	0.30
					
<b>KM5120-10-230</b>	Round Fluted Extended	mm	96.0	3.0	3.0
		in	3.78	0.12	0.12
<b>KM5120-10-240</b>	Round Fluted Extended	mm	96.0	4.0	4.0
		in	3.78	0.16	0.16
<b>KM5120-10-250</b>	Round Fluted Extended	mm	93.0	5.0	5.0
		in	3.66	0.20	0.20
<b>KM5130-10-230</b>	Round Fluted Extended	mm	149.0	3.0	3.0
		in	5.87	0.12	0.12
<b>KM5130-10-240</b>	Round Fluted Extended	mm	149.0	4.0	4.0
		in	5.87	0.16	0.16
<b>KM5130-10-250</b>	Round Fluted Extended	mm	146.0	5.0	5.0
		in	5.75	0.20	0.20

**TAVA Surgical® High Speed Replacement Burs**  
**Compatible with Stryker® TPS™, Medtronic® Midas Rex®**  
**Classic I & II and Linvatec® UltraPower™ Systems**



Order Number	Description	Unit	Overall Length	Head Diameter	Head Length
					
<b>KM-M340DU</b>	Round Diamond	mm	84.0	4.0	4.0
		in	3.31	0.16	0.16
					
<b>KM-U295</b>	Oval Aggressive Cut	mm	89.0	4.0	8.0
		in	3.50	0.16	0.31
<b>KM5110-15-40</b>	Oval Aggressive Cut	mm	72.0	4.0	9.5
		in	2.83	0.16	0.37
<b>KM5130-15-40</b>	Oval Aggressive Cut	mm	142.0	4.0	9.5
		in	5.59	0.16	0.37
					
<b>KM-M8U</b>	Match Head	mm	93.0	2.9	3.8
		in	3.66	0.11	0.15
<b>KM-AM8U</b>	Match Head	mm	147.0	2.9	3.8
		in	5.79	0.11	0.15
					
<b>KM-M52U</b>	Metal Cutter (carbide)	mm	95.0	3.2	17.0
	(Round End)	in	3.74	0.13	0.67

# The TAVA Surgical Value

TAVA Surgical has set the standard for excellence in offering surgical products compatible with original manufactures power systems.

- ✓ Low Cost / High Value – Manufacturer Direct
- ✓ Single Source Provider - Product Consolidation / Standardization
- ✓ GPO / IDN Contracts
- ✓ Products Manufactured in USA
- ✓ Product Cross Referencing on website
- ✓ Manufacturer Direct National Sales Force & Inside Account Managers
- ✓ Inquire about our Free 2nd Day Express Shipping.
- ✓ Customer Service - Product Assistance Help Desk  
(Extended Customer Service Hours 7:30PM EST)
- ✓ Quality Management System – ISO Certified



One Brasseler Boulevard • Savannah, Georgia 31419 • 912-921-7575  
800-569-6738 Ext. 7050 • 912-921-7578 (fax) • [www.TavaSurgical.com](http://www.TavaSurgical.com)

Linvatec® and UltraPower™ are trademarks of ConMed Corporation. Medtronic®, Midas Rex® and Legend® are registered trademarks of Medtronic, Inc.  
Stryker® and TPS™ are trademarks of Stryker Corporation.  
TAVA Surgical® is not affiliated with any of the above.