

Pulse Lavage

System & Accessories



Specifications ■

SPECIFICATIONS:

PULSE RATE

– Variable to 1,700 pulses per minute by trigger and pressure adjustment

FLOW RATE

– Variable to 0.9 liters per minute

OPERATING PRESSURE

– 75 to 110 psi (5-7 bar)

- Ergonomically designed and balanced, the L31 Pulse Lavage is an easy-to-use, pneumatic instrument that pumps fluids in distinct pulses at optimum pressure to remove blood and debris.

The L31 Pulse Lavage features an easy-to-load pump cartridge to facilitate easy operation.

A wide range of accessories including disposable tubing and nozzle sets designed for use during general purpose, hip and knee replacement surgeries are available.

Unique groove in the Pulse Lavage handle keeps the delivery tube neatly tucked in and out of the way.

Available with many popular hose connectors.



▲ Shown above with the knee arthroplasty set.


▼ Shown below with the 10.25"(25cm) femoral nozzle and suction attachment.



Attachment Sets

L21-030


General Purpose Set



- Pump Cartridge
- 3m of tubing
- Bag spike connector
- 10cm shower head

L21-031


Hip Arthroplasty Set



- Pump Cartridge
- 3m of tubing
- Bag spike connector
- 10cm shower head
- acetabulum nozzle
- 25cm femoral shaft nozzle

L21-032

Knee Arthroplasty Set



- Pump Cartridge
- 3m of tubing
- Bag spike connector
- 10cm angled shower head

Product Number Hose Connection

L31-100	Compatible to a Tava® Hose Connection
L31-200	Compatible to a Linvatec® & Hall® Hose Connection
L31-300	Compatible to a 3M® Hose Connection
L31-400	Compatible to a Stryker® Hose Connection
L31-500	Compatible to a Synthes® Hose Connection
L31-600	Compatible to a Aesculap® Hose Connection



Accessories

A90-000	Air Hose compatible to a Tava Surgical™/MicroAire® Connector (Length: 3m/10 feet)
A90-000-15	Air Hose compatible to a Tava Surgical™/MicroAire® Connector (Length: 5m/15 feet)
A90-107	Air Hose (UK MK7 Style) compatible to a Tava Surgical™/MicroAire® Connector (Length: 3m/10 feet)
A90-107L	Air Hose (UK MK7 Style) compatible to a Tava Surgical™/MicroAire® Connector (Length: 5m/15 feet)
A90-200	Air Hose compatible to a Linvatec® & Hall® Connector (Length: 3m/10 feet)
A90-001	Cleaning Plug/Cap compatible to a Tava Surgical™/MicroAire® Connector
A90-201	Cleaning Plug/Cap compatible to a Linvatec® & Hall® Connector
A90-501	Cleaning Plug/Cap compatible to a A.O. Synthes® Connector
A90-601	Cleaning Plug/Cap compatible to a Aesculap® Connector
M001	Instrument Lube


**Pulse Lavage
System & Accessories**



*Products not shown at
actual size.*

Order/Information _____ +1 912 921-7575 _ +1 912-925-8525

Instrument Accessories



L21-004

10cm Showerhead Nozzle with Locking Collar



L21-005

10cm Angled Showerhead Nozzle with Locking Collar



L31-042

Suction Tubing




L21-007

22.5cm Femoral Shaft Nozzle with Locking Collar



L31-040

Femoral Suction Attachment




L31-046

Suction Tubing and Femoral Suction Attachment



L21-010

Splash Shield, 15cm Round



L21-013

Splash Shield, 12.5cm x 20cm Rectangular



L21-017

Cone (Cup) Splash Shield Ø14cm



L21-020

Cone (Cup) Splash Shield Ø14cm



TAVA SURGICAL™
SURGICAL POWER & ACCESSORIES

One Brasseler Boulevard ▪ Savannah, Georgia 31419 ▪ +1 912-921-7575
+1 912-925-8525 ▪ +1 912-921-7578 (fax) ▪ www.tavasurgical.com

Order/Information

+1 912-921-7575 or +1 912-925-8525

Stryker® is a registered trademark of Stryker Corporation. Linvatec® and Hall® are registered trademarks of Linvatec Corporation. 3M® is a registered trademark of 3M Corporation. MicroAire® is a registered trademark of MicroAireSurgical Instruments, Inc. Synthes® is a registered trademark of Synthes (U.S.A.). Aesculap® is a registered trademark of Aesculap Instruments Corporation. *Tava Surgical™ is not affiliated with any of the above. Tava Surgical™ makes no claim to copyrights or trademarks which are the property of these companies.*

All rights of distribution, also by photocopy, reprint or storing and recovery in any kind of data processing are reserved and require our written approval. Color and products subject to alterations. Printed in USA.

