

Trepan Motor



The Tava Surgical™ Trepan motor, designed to accept cranial perforators with Hudson fittings, provides the surgeon an efficient, precise instrument with a motor rotation speed at 750 rpm – 850 rpm.

SPECIFICATIONS	Trepan Motor
Power	Compressed sterile air or nitrogen
Pressure Supply	120 psi - 130 psi
Motor Rotation Speed	850 rpm
Motor Length	171 mm
Motor Diameter	38 mm
Sterilization	Complete system can be hand wiped with a damp cloth and autoclaved (132°C/270°F) NOTE: Do Not Immerse!



Order/Information

912-921-7575 or 912-925-8525



TAVA SURGICAL™
SURGICAL POWER & ACCESSORIES

One Brasseler Boulevard • Savannah, Georgia 31419

912-921-7575 • 912-925-8525 • 912-921-7578 (fax)

Expect Quality. **Demand Tava Surgical™**