



TAVA SURGICAL™
SURGICAL POWER & ACCESSORIES

One Brasseler Boulevard ▪ Savannah, Georgia 31419 ▪ +1 912-921-7575
+1 912-925-8525 ▪ +1 912-921-7578 (fax) ▪ www.tavasurgical.com

All rights of distribution, also by photocopy, reprint or storing and recovery in any kind of data processing
are reserved and require our written approval.

Color and products subject to alterations.


Printed in USA.

Order/Information
+1 912-921-7575 or +1 912-925-8525



Tava Surgical™ – Surgical Power & Accessories

XK-PRO 100™ High Speed Cutters



Expect Quality. Demand Tava Surgical™.



TAVA SURGICAL™
SURGICAL POWER & ACCESSORIES



Tava Surgical™

XK-PRO 100™
High Speed
Cutters

CUTTER COMPATIBILITY


Straight and Angled Attachments / Cutters

Each group of nosepiece attachments (Short 50 mm, Medium 70 mm, Long 90 mm, Extra Long 110 mm) are designed to operate with a wide range of dissecting cutters. Straight and Angled attachments of 50, 70, 90 and 110 mm lengths must be used with 50, 70, 90, 110 mm length cutters respectively.

Example:

Cutter Item number S30RDX indicates the following:

- S = 50 mm attachment
- 30 = 3.0 mm head diameter
- R = Round head shape
- D = Diamond head type
- X = Extra Coarse



Note: Each length cutter must only be used with the corresponding matching nosepiece attachment.


Craniotome Attachments / Cutters


Each Craniotome attachment must only be used with its corresponding size cutter.

Size	Craniotome Attachment	Craniotome Cutter
Pediatric	XKP-C12	CF12S
Adult	XKP-C16	CF16S



Tava Surgical™
XK-PRO 100™
High Speed
Cutters

Description	Attachments/Catalog Number				Head Size mm	No. of Flutes or Grit
	S (50mm)	M (70mm)	L (90mm)	X (110mm)		
 Round Fluted	S10RF	M10RF	-	-	Ø 1.0	6
	S20RF	M20RF	-	-	Ø 2.0	6
	S30RF	M30RF	L30RF	X30RF	Ø 3.0	8
	S40RF	M40RF	L40RF	X40RF	Ø 4.0	8
	S50RF	M50RF	L50RF	X50RF	Ø 5.0	8
	S60RF	M60RF	L60RF	X60RF	Ø 6.0	10
	S70RF	M70RF	-	-	Ø 7.0	12

 Round Diamond	S06RDF	M06RDF	-	-	Ø 0.6	Fine
	S10RDF	M10RDF	L10RDF	X10RDF	Ø 1.0	Fine
	S20RDC	M20RDC	L20RDC	X20RDC	Ø 2.0	Coarse
	S30RDX	M30RDX	L30RDX	X30RDX	Ø 3.0	X Coarse
	S40RDX	M40RDX	-	-	Ø 4.0	X Coarse
	S50RDX	M50RDX	-	-	Ø 5.0	X Coarse
	S60RDX	M60RDX	-	-	Ø 6.0	X Coarse
	S70RDX	M70RDX	-	-	Ø 7.0	X Coarse

 Acorn	S50AF	M50AF	-	-	Ø 5.0 x 7.0	8
	S60AF	M60AF	-	-	Ø 6.0 x 7.5	8

 Match Head	-	M30MF	-	X30MF	Ø 3.0 x 3.8	2
---	---	-------	---	-------	-------------	---

 Pin Cutter	-	M10PC	-	-	Ø 1.0 x 8.0	2
---	---	-------	---	---	-------------	---

Description	Attachments/Catalog Number		Head Size mm	No. of Flutes or Grit
	Attachment	Cutter		
12mm	XKP-C12	CF12S	Ø 1.8 x 12.0	2
16mm	XKP-C16	CF16S	Ø 1.8 x 16.0	2